Unimatrix R

Hard ◦ Soft

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INSTRUCTIONS FOR USE (EN)

The Unimatrix R System is indicated for Class II restorations and composed by preformed sectional matrices, available in 3 sizes, including rings for its attachment and silicone tips. The main advantages of this system are the easiness in the reconstruction of contact points and natural tooth anatomy, and comfort to the patient. The ends of the ring are endowed with triangular and anatomical silicone tips, optimizing the adaptation of the matrix. The Unimatrix R clamps are available in two models: Hard (blue) indicated for cases which require the production of contact points with greater dental separation and pressure and the Soft clamp (pink) indicated for clinical cases which require less dental separation and pressure.

Use: Isolate the operative field with a rubber dam or other suitable isolation technique. Select, among the 3 matrix sizes, the item which is best adapted to the cavity to be restored. Install the matrix in the proximal surface to be restored. When it concerns a MOD cavity, it is advisable to restore one surface at a time. Install the wedge. TDV elastic wedges improve the adaptation of the matrix. In this case, silicone tips are not used and the ring is positioned in the elastic wedge piercings. Install the Soft (pink) or Hard (blue) Unimatrix R ring as per the clinical case. For opening and inserting the ring, TDV Unimatrix pliers or Palmer pinchers can be used. The use or not of the ring and silicone tips and the means of their installation depend upon personal preferences and the clinical situation. Having concluded the restoration, remove first the fastening clamp, then the wedge and lastly the matrix. To obtain consistent contact points, it is advisable to separate the adjacent ones by means of wedges or fastening ring, from the anesthetic until the installation of the matrix.

Physical principle of functioning: The matrix acts as a physical barrier for the restorative material. Its format allows the anatomical reconstruction of the tooth. The ring is made of a metallic alloy which has memory effect and sufficient force to fasten the matrix and separate the teeth.

Warnings, precautions and restrictions: Professional use. The matrices are of single use. Their reuse can cause cross-contamination. When there is proven allergy to any one of its components, the product must not be used. Excessive expansion with the pliers can lead to deformation of the ring. During the insertion of the ring, maintain the finger pressed on it until the removal of the pliers, ensuring the adjustment and avoiding sliding. Care must be taken during the insertion of the ring, due to its flexibility. Ensure that the product is totally fitted in the pliers. Avoid leaving the ring open in the pinchers for a long period of time. It is advisable to protect the papilla by any method to avoid accidents. It is advisable that the use of the product is always accompanied by absolute isolation of the operative field, as objects can be swallowed or breathed in by accident. The ring must be used with the silicone protectors for an appropriate adjustment of the matrix. The ring and the silicone tips must be sterilized separately.

Sterilization: The ring and silicone tips are sterilizable in autoclave until 134 °C - 18 min. Before sterilization, execute the cleaning with the aid of water and enzymatic detergent. The sterilization can cause slight change in the coloring of the protectors and reinforcement, but does not affect their functionality. Composition: Matrizes: stainless steel. Ring: nickel-plated carbon steel. Ring reinforcement: thermoplastic. Tips: silicone.

Ref.	Code	
04126H	0003020081	Unimatrix R ring Hard – 1 ring + 4 triangular silicone tips
04124R	0003020085	Unimatrix R – 50 matrizes (S/M/L) + 2 rings + 8 triangular silicone tips
04122R	0003020086	Unimatrix R – 25 matrizes (S/M/L) + 1 ring + 4 triangular silicone tips
04120R	0003020087	Unimatrix R – 100 matrizes (S/M/L) + 2 rings + 8 triangular silicone tips
04126S	0003020102	Unimatrix R ring Soft – 1 ring + 4 triangular silicone tips
04159T	0003020103	Triangular silicone tips for Unimatrix rings – 20 units
04160T	0003020104	Triangular silicone tips for Unimatrix rings – 50 units

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(EN) Preformed sectional matrices system for Class II restorations

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Obelis S.A