

Product catalogue



FSC



TePe is certified according to the international quality standard ISO 9001 and the environmental standard ISO 14001.

TePe Munhygienprodukter AB
 Bronsåldersgatan 5 | 213 76 Malmö, Sweden
 Tel +46 (0)40 670 11 00 | Fax +46 (0)40 670 11 11
 info@tepe.com | www.tepe.com

All trademarks herein are the property of TePe Munhygienprodukter AB.



© TePe Munhygienprodukter AB 2025. 1321BET

Table of contents

Tiny things for big smiles	4
Prevention	5
Sustainability	6
Interdental cleaning	7
Toothbrush	15
Speciality brushes	23
Healthy smiles	27
Education	33
Customised products	35
Marketing material	36



Some products have symbols indicating their suggested use and available packaging types. Here is a reference guide for all symbols:

Areas of use:	Packaging:	
Basic oral care for those with healthy mouths	Flow pack	
Suitable for children and/or juniors	Blister	Bag
Suitable for those with braces	Carton pack	Box
Suitable for those with implants	Tube	Bottle
Suitable for those with gingivitis/periodontitis		

Reference:

- Greenhouse gas emissions are calculated in accordance with the GHG Protocol and ISO 14040 and ISO 14044.
- TePe sources 100% raw materials for the toothbrush- and interdental brush handles, as well as the flosser, tongue cleaner, travel case, cap and extra grip, from either renewable raw materials or certified renewable raw materials via the mass balance approach, the latter independently certified by the International Sustainability and Carbon Certificate PLUS (ISCC PLUS).
- Furness S, et al. Interventions for the management of dry mouth: topical therapies. Cochrane Database Syst Rev. 2011 Dec 7;(12):CD008934.
- Key findings within the group saying dry mouth is affecting their daily life (13%). TePe – Mouth Dryness in UK, web panel questionnaire study with 1000 respondents, Transvector, April 2023.


Tiny things for big smiles

It's often the tiny things in life that make it meaningful, like starting the day by sharing a smile with a friend, colleague or a stranger. The same is true for our oral health; a few minutes of toothbrushing, flossing, and in-between brushing can greatly affect our overall health and wellbeing.

It's a tiny effort, but a huge deal.

Taking care of your mouth means taking care of your health. That's why we continue innovating, educating, and providing a broad range of products for better oral health.

We started in 1965, and have continued to find clever solutions ever since. Look through this catalogue and discover the tiny things TePe offers to make people lead healthier lives with healthy smiles.



Prevention starts with daily care

You have probably heard us talk about prevention many times. That's because it's really important to us. Caring for your teeth and gums is more than just a routine to make your mouth feel fresh; it's essential for preventing common issues like cavities and gum disease. In other words, brushing, flossing, and regular

use of mouthwash – can save you from more extensive and expensive treatments in the future. Our products are designed to make daily prevention simple, effective, and enjoyable. Investing a little bit in our oral health today, can contribute to a healthier tomorrow.

1.



Education about prevention

TePe raises awareness about the positive impact of prevention with the aim to...

- Reduce development of oral diseases
- Reduce reoccurrence of oral diseases
- Improve or maintain level of oral diseases

2.

Increased knowledge improves health

Preventive care leads to healthier mouths which supports general health through...

- Functional wellbeing
- Psychological and social wellbeing
- The absence of discomfort or pain



3.



Improved health is good for the planet

Good general and oral health has a positive impact on the environment since it can lead to...

- Fewer resources used
- Less waste generated
- Reduced emissions

Sustainability is the healthy approach

Taking care of your teeth is great for your smile. Did you know that it can also contribute to creating a healthier planet? Healthy teeth often mean a lower environmental impact, as it can reduce the need for treatments and the emissions associated with dental repairs.

What about our products and production? In Sweden, all of our products are made using 100% renewable energy, boosted by our own solar park. We actively work to reduce greenhouse gas emissions and are phasing out virgin fossil plastics in favour of renewable options.^{1,2} Sustainability is a journey where improvements are made step by step, where every action counts – both big and small. That is what we call a healthy approach. Take a look at some of our milestones.

This is the beginning of our sustainability journey. TePe is founded with the launch of the triangular wooden dental stick, which is made of Swedish wood and developed in collaboration with the Faculty of Odontology in Malmö, Sweden.

1965



TePe starts using hydro power. The production site is now powered by 100% renewable energy.

1993

Innovation of preventive care: the TePe Interdental brush is launched.

2016



2019




Starts to phase out virgin fossil material. Introducing a Life Cycle Assessment for self-produced products and packaging.



A solar panel park is installed to power our production.

2018



TePe's cross-functional sustainability team is established.

Honeybees find a new home on our roof.

2022

All the handles for the TePe tooth-brushes and interdental brushes are from now on made with renewable or mass balance materials.




2021

Linking 7 of the Global Goals to our strategy.



2030

Contribute to 25 million healthy smiles.



2023

Launching a tooth-brush with wooden handle and replaceable brush head and a dental floss made from recycled plastic bottles.



2024

Greenhouse gas mapping of scope 1, 2 and 3.



2025

100% of the packaging should now include renewable, recycled or FSC-certified materials.





Interdental cleaning

A special focus on interdental cleaning will keep those gaps gleaming and is essential for maintaining a healthy smile. We offer a range of interdental products designed to suit every individual and every type of gap.

Recommend the correct interdental cleaning product

When there are enough spaces between teeth, an interdental brush is the ideal tool. It can effectively reach and clean the nooks and crannies of the teeth, including the concave surfaces that floss may miss.

In narrow spaces where the gums completely fill the gap between the teeth, using floss or a toothpick is often sufficient to keep the area clean.




TePe® Interdental Brush Original

According to studies, interdental brushes are the most effective tool for removing plaque between teeth. Our brushes are developed together with dental experts, to ensure the highest functionality, comfort and quality. We offer 9 colour-coded ISO standard sizes, so there is a size for every type of gap, both big and small. Matching the gaps with the right colour makes it easy to remember what size you need. All the handles of our interdental brushes are made with renewable materials through the mass balance approach.²

The bristles cover the entire wire for safe and gentle cleaning.

Plastic-coated wire and rounded tip for safety and comfort.



The smaller brush sizes (0-3) have a flexible neck for increased durability and access.



The larger brush sizes have a thicker, bendable wire.



Ergonomic handle for a stable grip.



It is user-friendly thanks to the colour-coded sizes.



The plastic cap can be used for hygienic protection or to extend the handle.





TePe® Interdental Brush Original is also available in a mixed pack.





TePe® Interdental Brush Original



Colour: Pink
ISO size: 0
PHD size: 0.6 mm
Wire size: 0.4 mm




Colour: Orange
ISO size: 1
PHD size: 0.7 mm
Wire size: 0.45 mm





Colour: Red
ISO size: 2
PHD size: 0.9 mm
Wire size: 0.5 mm





Colour: Blue
ISO size: 3
PHD size: 1.1 mm
Wire size: 0.6 mm





Colour: Yellow
ISO size: 4
PHD size: 1.3 mm
Wire size: 0.7 mm




Colour: Green
ISO size: 5
PHD size: 1.6 mm
Wire size: 0.8 mm




Colour: Purple
ISO size: 6
PHD size: 2.0 mm
Wire size: 1.1 mm



Colour: Grey
ISO size: 7
PHD size: 2.6 mm
Wire size: 1.3 mm



Colour: Black
ISO size: 8
PHD size: 3.4 mm
Wire size: 1.5 mm




TePe® Interdental Brush Extra Soft


- Extra soft bristles for gentle and safe cleaning.
- Six sizes, all with plastic-coated wire.
- For patients with sensitive teeth, sore gingiva, dry mouth, or after oral surgery*.




Colour: Orange
ISO size: 1
PHD size: 0.7 mm
Wire size: 0.45 mm




Colour: Blue
ISO size: 3
PHD size: 1.1 mm
Wire size: 0.6 mm




Colour: Green
ISO size: 5
PHD size: 1.6 mm
Wire size: 0.8 mm




Colour: Red
ISO size: 2
PHD size: 0.9 mm
Wire size: 0.5 mm




Colour: Yellow
ISO size: 4
PHD size: 1.3 mm
Wire size: 0.7 mm




Colour: Purple
ISO size: 6
PHD size: 2.0 mm
Wire size: 1.1 mm



Colour: Pink
ISO size: 0
PHD size: 0.6 mm
Wire size: 0.4 mm




Colour: Red
ISO size: 2
PHD size: 0.9 mm
Wire size: 0.5 mm




Colour: Yellow
ISO size: 4
PHD size: 1.3 mm
Wire size: 0.7 mm



Colour: Orange
ISO size: 1
PHD size: 0.7 mm
Wire size: 0.45 mm



Colour: Blue
ISO size: 3
PHD size: 1.1 mm
Wire size: 0.6 mm



Colour: Green
ISO size: 5
PHD size: 1.6 mm
Wire size: 0.8 mm



- Plastic-coated wire for safe use.
- Long, flat handle for a stable grip.
- Angled brush head for easier access, especially between the posterior teeth.



TePe Angle™ is also available in a mixed pack.

* A dental professional should always be consulted regarding brushing after oral surgery.



TePe Interdental Brush Extra Soft is also available in a mixed pack.



TePe EasyPick™ – easy cleaning in between

This soft silicone toothpick efficiently removes plaque and pesky food debris to keep your teeth and gums healthy and clean. It's flexible and comfortable to use, even if you have braces, implants, or crowns.

Unique design and silicone coating ensure efficient and comfortable cleaning.

Rounded tip for comfortable use.


The unique lamellae provide a wide contact area with the surfaces between the teeth.

Overlapping conical sizes to fit narrow as well as wide spaces between teeth.

A handy pocket case for TePe EasyPick™ is included for easy cleaning on-the-go.

TePe EasyPick™

This product is a good choice for healthy mouths, in the case of gingivitis, and as an on-the-go alternative for everyone.



Colour: Orange

Size: Extra Small/Small



Colour: Blue

Size: Medium/Large




Colour: Purple

Size: Extra Large



TePe Mini Flosser™

- Pre-loaded dental floss holder for easy and convenient use in narrow spaces.
- Strong, slim and shred-resistant floss.
- Bite plane for safe and gentle insertion.
- Handle made with renewable materials through the mass balance approach.²




Insert the floss by gently biting onto the bite plane.

TePe® Gingival Gel

Size: 20 ml

- Contains chlorhexidine (0.2%) to fight bacteria, and fluoride (0.32% NaF) to strengthen the tooth enamel and protect sensitive root surfaces.
- Non-abrasive, without alcohol and colourants.
- Fresh mint flavour and smooth formula.
- Provides extra caries protection.




Use the cap as a bowl to apply gel.



TePe® Dental Floss

Length: 5 m / 40 m

- Floss is recommended for tight interdental spaces, where an interdental brush will not fit.
- Expanding floss that suits different gaps as it adapts to the dental space.
- Slides easily between teeth due to its vegetable wax and avocado oil coating.
- Fresh mint flavour.



Gently move the floss along each tooth surface.



TePe® Travel Case


- Practical travel case for TePe® Interdental Brush, TePe Mini Flosser™ or TePe EasyPick™.
- Air holes to allow bristles to dry.
- For hygienic storage during travel.
- Made with renewable materials through the mass balance approach.²



TePe® Bridge & Implant Floss

Amount: 30 pcs

- Wide and spongy floss for efficient cleaning.
- The short length and sturdy ends ease insertion and use.
- Invaluable for cleaning of implants, bridges and braces.
- Daily use contributes to healthy gums.





Toothbrushes

What all bristles have been missing. Now everybody can put TePe toothpaste on their TePe brushes. The TePe Daily & Pure range offers mild fluoride toothpaste for every family member and for people with a sensitive or dry mouth.

A simple 3-step routine for a healthier smile

Always start by cleaning between your teeth to remove plaque. Using either an interdental brush or a dental floss is actually the most effective way to remove plaque between the teeth.



The tongue cleaner is the magic wand against bad breath! The TePe Tongue Cleaner has three blades for a triple cleaning effect. It's easy to use, removes bacteria with just a few strokes and prevents bad breath.




Start by brushing the insides of your teeth to make sure they get as much love as the outsides. Brush twice a day for two minutes for optimal plaque removal and fluoride distribution.





TePe Supreme™ – supreme plaque removal with a gentle feel

Our innovative toothbrush that offers a longer reach, thanks to the two-level bristle design. The longer bristles clean further between the teeth and along the gum line, while the shorter bristles gently clean the rest of the tooth surfaces. Developed in collaboration with dental experts. And by the way, all the handles of our toothbrushes are made with renewable materials through the mass balance approach.²



Unique two-level bristle design for excellent reach.



Easy-grip handle with thumb pad for extra comfort.

Can be angled in hot water for better access.


End-rounded and tapered filaments.



TePe Supreme™

Size: Regular
Bristles: Soft


- Unique two-level bristle design for excellent reach.
- Easy-grip handle with thumb pad for extra comfort.



TePe Supreme™ Compact

Size: Compact
Bristles: Soft


- Unique two-level bristle design for excellent reach and easy-grip handle with thumb pad.
- Suitable for adolescents or for adults who prefer a smaller brush head.



TePe Nova™

Size: Regular
Bristles: Extra soft (pink), soft (yellow) & medium (blue)

- Long brush tip and end-rounded bristles for effective and gentle cleaning.
- Easy-grip handle with thumb pad for extra comfort.



The 4x2 rule is worth sharing



Use about 2 cm of fluoride toothpaste for adults.

(Children 0–2 years: a rice grain-sized amount.
3–6 years: a pea-sized amount)



Brush your teeth morning and evening.



Brush for at least two minutes.





Avoid eating or drinking for two hours to get the most benefit from the fluoride.

TePe Select™ Regular

Size: Regular
Bristles: Extra soft, soft & medium

- Tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- User-friendly handle design.



TePe Select™ Compact

Size: Compact
Bristles: Extra soft, soft & medium


- Small, tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- Suitable for children and adults who prefer a smaller brush head.



TePe Mini™ Bolibompa

Size: Mini
Bristles: Extra soft


- Extra small, tapered brush head and end-rounded bristles for gentle and effective cleaning.
- Recommended from the first tooth.



TePe Mini™

Size: Mini
Bristles: Extra soft & soft


- Extra small, tapered brush head and end-rounded bristles for gentle and effective cleaning.
- Recommended from the first tooth.
- Available with or without charming prints created by Swedish Isa Form.



TePe Colour™

Size: Regular
Bristles: Soft

- Tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- User-friendly handle design.



TePe Colour™ Compact

Size: Compact
Bristles: Extra soft


- Small, tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.
- Suitable for children or for adults who prefer a smaller brush head.



TePe Kids™

Size: Compact
Bristles: Extra soft & soft

- Toothbrush with colourful print, suitable from the age of 6 years.
- Small, tapered brush head for easy access.
- End-rounded bristles for gentle and effective cleaning.



A good routine for healthy teeth in the future

Kids need help with brushing their teeth

Until they're 10-12 years or so. Why not make it fun by adding a song or a story?



Create a habit

Brush the inside, outside and biting surfaces in the same order to make sure no side is forgotten. Remember the gum line where bacteria can easily build up and cause inflammation.






TePe Choice™ – brush away in a natural way

TePe Choice™ is an innovative and user-friendly toothbrush with a reusable handle made from Swedish beechwood, and three replaceable brush heads.

Size: Regular
Bristles: Soft


- Tapered brush head for easy access.
- Plant-based coating on the handle to minimise moisture absorption.
- Soft, end-rounded bristles for gentle and effective cleaning.
- Brush head made with renewable materials through the mass balance approach.²



TePe Classic™

Size: Regular
Bristles: Soft

- User-friendly handle design.
- End-rounded bristles for gentle and effective cleaning.
- For those who prefer a bigger brush head.



TePe PlaqSearch™

Size: 10 tablets/per pack

- A disclosing product for quick and simple self-check identification of plaque.
- Two-tone colour indicators.
- Teaching aid for better oral health.




Differentiates new (red) and old (blue) plaque.



TePe® Extra Grip

- Provides a comfortable, stable grip which helps patients with reduced manual dexterity.
- Fits the Select toothbrush and the specialty brushes.
- Easy to insert and remove by tapping the toothbrush handle on a hard surface.
- Weighs only 30 grams.
- Made with renewable materials through the mass balance approach.²



Insert/remove the brush by tapping the handle against a table.




TePe® Brush Head Cap

- Small, practical travel cap.
- Air-holes ensure quick drying of the brush.
- Suits all TePe toothbrushes and special brushes (except the Denture and Classic specialty Brush).
- Made with renewable materials through the mass balance approach.²



Speciality brushes


Wisdom teeth, dental implants, crowns/bridges or braces – no matter your needs, TePe can help make cleaning easier. TePe offers a wide range of high-quality cleaning devices with excellent functionality for those with specific oral hygiene needs.



TePe Special Care™

Size: Regular (red) & Compact (blue)
Bristles: Ultra soft

- For patients with extremely sore and delicate tissues, e.g. as a result of trauma, radiotherapy treatment or after oral surgery.
- 12,000 (Regular) and 9,000 (Compact) ultra soft bristles for delicate cleaning.
- First brush stage after oral surgery.*



* A dental professional should always be consulted regarding brushing after oral surgery.



TePe Gentle Care™

Size: Regular
Bristles: Super soft

- For patients with dry mouth, sensitive root surfaces or sensitive oral tissues.
- Super soft filaments for gentle cleaning.
- Second brush stage after oral surgery.*



TePe® Denture Care

Size: Regular

- Effective and safe cleaning of removable dentures.
- Brush head can be bent backwards for improved grip and accessibility.
- Designed to resemble a regular toothbrush.



Long, strong filaments for efficient cleaning of dentures.



TePe® Tongue Cleaner

- Three blades for a triple cleaning effect.
- Helps prevent bad breath (halitosis).
- Removes bacteria in a few strokes.
- Made with renewable materials through the mass balance approach.²



TePe Compact Tuft™

TePe Compact Tuft is an angled brush designed for precision. It features a small, dome-shaped brush with tightly packed short filaments. With its stable handle and dense, firm tuft, this brush makes it easier to clean even the hardest-to-reach areas. And by the way, all the handles of our specialty brushes are made with renewable materials through the mass balance approach.²

Dome-shaped brush head with dense, rounded tuft.

Tightly packed short bristles that give a firm feel.




TePe Compact Tuft™

- Optimal for precision cleaning.
- Dome-shaped brush head with dense, rounded tuft.
- Ideal for cleaning along the gum line, on the occlusal surfaces of emerging teeth, or impacted wisdom teeth.



TePe Universal Care™

- A unique angle provides easy access to the lingual and palatal surfaces of teeth and implants.
- Makes it easier to clean narrow and hard-to-reach areas.
- End-rounded bristles for gentle cleaning.





The angled neck makes it easier to clean the lingual surfaces.



TePe Interspace™

Bristles: Soft (yellow) & medium (blue)

- Ideal for cleaning furcations, implants, fixed braces or behind the last tooth.
- Brush tip can be attached to either side of the handle.
- Comes with 12 replaceable brush tips.



TePe® Implant Orthodontic Brush

- Extra slim brush head for improved access.
- End-rounded bristles for gentle and effective cleaning.
- Optimal for cleaning buccal implant surfaces or orthodontic appliances.
- Long and thin neck for better access.



Effective cleaning around brackets and wire.



TePe® Implant Care Kit*

A collection of carefully selected maintenance products designed for effective plaque removal in and around implants, making the cleaning simple.

- Products included:
- TePe Compact Tuft™
 - TePe Implant Orthodontic™ Brush
 - TePe Universal Care™
 - TePe® Bridge & Implant Floss (30 pcs)
 - TePe® Interdental Brushes Extra Soft mixed pack
 - Patient Information Leaflet



TePe® Orthodontic Kit*

The kit features carefully selected self-care products and readily accessible tips and advice for orthodontic patients.

- Products included:
- TePe Supreme™ Compact
 - TePe Compact Tuft™
 - TePe® Interdental Brush, ISO size 4 and 5
 - TePe Angle™ Interdental brush, ISO size 4
 - TePe EasyPick™, XS/S and M/L
 - TePe® Orthodontic wax
 - Patient Information Leaflet



*Currently available in selected markets.





Healthy smiles for all mouths

TePe offers a wide selection of products suitable for all stages of life – from the first teeth to special needs like oral care for braces and implants. You will find our product recommendations throughout the following pages.



Paedodontic care

To keep the smile bright and healthy for a lifetime, we suggest brushing as soon as that first tooth pops up. It is often recommended to start with interdental cleaning during the teenage years.



TePe Supreme™ Compact

TePe Compact Tuft™

TePe Select™ Compact

TePe Mini™

TePe PlaqSearch™

This section displays paedodontic products. At the top, two circular diagrams illustrate the brush heads: 'TePe Supreme™ Compact' (orange) and 'TePe Compact Tuft™' (blue). Below these are a yellow brush, an orange brush, a blue brush, and a green brush. In the center is a box of 'TePe PlaqSearch™' interdental brushes with three purple brushes next to it. At the bottom, two circular diagrams show 'TePe Select™ Compact' (orange) and 'TePe Mini™' (blue) brushes being used on teeth. The 'TePe PlaqSearch™' label is positioned below the central box.

The foundation of maintaining healthy gums is daily interdental cleaning and toothbrushing. Professionals play a key role in motivating patients and recommending the right products.

Daily care



TePee® Tongue cleaner

TePe Mini Flosser™

TePe Supreme™

TePe® Dental Floss

TePe Angle™

TePe EasyPick™

TePe® Gingival Gel

TePe® Interdental Brushes

This section displays adult dental care products. At the top, four circular diagrams illustrate the products: 'TePee® Tongue cleaner' (green), 'TePe Mini Flosser™' (green), 'TePe Supreme™' (orange), and 'TePe® Dental Floss' (white). Below these are a green tongue cleaner, a green flosser, a yellow brush, and a roll of white dental floss. In the center is a green bottle of 'TePe® Gingival Gel'. At the bottom, four circular diagrams show 'TePe Angle™' (green), 'TePe EasyPick™' (orange), 'TePe® Gingival Gel' (green), and 'TePe® Interdental Brushes' (blue) being used on teeth. The 'TePe® Gingival Gel' label is positioned below the central bottle.



Implant care

Getting the right advice on suitable products and cleaning techniques can make a big difference in keeping your implants healthy. Daily cleaning helps to prevent peri-implant diseases.



TePe® Interdental Brushes



TePe Compact Tuft™

TePe® Implant Orthodontic




TePe Universal Care™

TePe® Bridge & Implant Floss




TePe EasyPick™

TePe Universal Care™




TePe® Interdental Brushes

TePe Supreme™ Compact



TePe® Implant Orthodontic

TePe Interspace™



TePe Compact Tuft™

TePe PlaqSearch™



TePe EasyPick™

Orthodontic care



Fixed orthodontic appliances, like braces, can quickly accumulate plaque around brackets, bands, and wires. To keep the smile healthy and bright, it's important to establish a good self-care routine that helps prevent gingivitis and caries.



Our global Clinical and Education Team:

RDH Lina Gassner Kanter

Clinical and Education Advisor

Dr Michaela von Geijer

Clinical and Education Advisor

Dr Sanjay Haryana

Clinical and Education Advisor

We provide clinicians with knowledge, tips, and tricks to achieve better compliance and treatment results.



Want a bigger bite of science?

Sign up for TePe Clinical Update

Four times a year, we will provide dental professionals and students with a bundle of inspiration to apply in clinical practice. We offer our newsletter Simply Science, a Webinar and Tips and Tricks.

Simply Science is a valuable resource keeping you updated with the latest scientific research in oral health. It is a brief summary of a recent scientific article, along with key insights that you can directly apply to your clinical practice. Our Clinical and Education Team will focus on sharing topics related to

TePe's core mission and the work of dental professionals - preventing oral disease and promoting healthy smiles worldwide.

We will continue with our highly appreciated Webinars. Rest assured, as long as you register, you can access the recording even if you don't attend the live session. We don't want anyone to miss out on valuable knowledge.

[Explore our knowledge platform at tepe.com](https://tepe.com)



Articles



Lectures



Webinars



Newsletter



Events



Sharing materials




Scan the QR-code to sign up



Practice Box

- Practical storage of TePe's interdenal brushes in the dental practice.
- 8 removable trays for IDB Original and Extra Soft, and 4 trays for TePe Angle. Each tray holds 15 brushes.
- The box includes a helpful chart that displays various brush sizes.



Customised products

Put the clinic's logo and a tagline on a toothbrush and spread awareness about the services, while promoting better oral health. TePe's imprinting service is an excellent marketing opportunity for dental practices.



Logotypes, images and messages can be printed on our Select, Select Compact, Colour, Colour Compact and Mini toothbrushes. Choose between a variety of different handle colours and print colours for amazing designs.



TePe® Dental Model

- For patient demonstration and instruction.
- Shows crowding, various interdental spaces, furcation, an erupting wisdom tooth and tooth loss.
- Can be marked with a non-permanent pen for instructional use.
- Suitable for the product cases featured below.



x1 Delivered in a brown hard paper box.

TePe® Product Demonstration Case

Prepacked with the majority of our products to show the variation of the product range. Includes extra space for marketing material and samples. Closes with a zipper. Carry it with ease via the shoulder strap or handle. Suited for sales representatives.



TePe® Product Presentation Case


Case featuring the major part of the TePe range. Suitable for dental clinics, for example.



Marketing & display material

TePe offers a range of information and demonstration materials to aid oral hygiene advice. Our product racks are perfect for display and storage.

Before ordering, please reach out to your TePe contact person, to make sure the design and product is available on your market.



Floor Stand, Metal

Item no: 931 054
Floor stand on wheels.
21 hangers, 1 shelf.
2 leaflet holders.



Floor Stand, Carton

Item no: 932 092
Floor stand.
4 shelf.



Mini Flosser Dispenser Box

Item no: 622 570 (Mini flosser, 120 pcs)
642 571 (Dental floss, 80 pcs)
Dispenser box for Mini flosser / Dental floss samples.



Medium Table Display

Item no: 932 033
Table display in cardboard.
8 hangers.



Wall Display, Metal

Item no: 931 053
Wall display
12 hangers. 1 shelf.
(Plastic containers available separately, item no. 931 038).



Table Display, Metal

Item no: 931 055
Table display with rotating function.
8 hangers.



Table Display, Carton

Item no: 932 091
Fits 2 SRPs.
(Toothpaste holder available separately, item no. 932 094).

Patient leaflet Display

Item no: 931 050
Display with 3 pockets for patient folders.
Clear acrylic material.



Mini Display

Item no: 931 031
For 6 toothbrushes.
9 interdental brushes.
Clear acrylic material.



Paper Bag

Item no: 256 INT



Toothbrush Display

Item no: 931 029
For 13 toothbrushes (Select range).
Clear acrylic material.



Interdental Brush Display

Item no: 931 035
For 15 interdental brushes.
Clear acrylic material.



EasyPick Display

Item no: 221INT
Cardboard table display.
For 3 SRPs, EasyPick.
(shelf ready packaging)



Sample Box

Item no: 152
Filled with a variety of product samples and information.

*Please note the visual design might change.